



Scala di lettura del cilindro graduato da 250 ml (millilitri).


Leggere le frazioni di divisione, i decimi di divisione.



Quanto vale questo livello ?
e poi gli altri scendendo




$$250\text{ml} - 1,5\text{div} = 250\text{ml} - 1,5 \cdot 2\text{ml} = 250\text{ml} - 3\text{ml} = 247\text{ml}$$



$$250\text{ml} - 1,5\text{div} = 250\text{ml} - 1,5 * 2\text{ml} = 250\text{ml} - 3\text{ml} = 247\text{ml}$$


$$230\text{ml} - 0,4\text{div} = 230\text{ml} - 0,4 * 2\text{ml} = 230\text{ml} - 0,8\text{ml} = 229,2\text{ml}$$



$$250\text{ml} - 1,5\text{div} = 250\text{ml} - 1,5 * 2\text{ml} = 250\text{ml} - 3\text{ml} = 247\text{ml}$$

$$230\text{ml} - 0,4\text{div} = 230\text{ml} - 0,4 * 2\text{ml} = 230\text{ml} - 0,8\text{ml} = 229,2\text{ml}$$

$$210\text{ml} + 0,1\text{div} = 210\text{ml} - 0,1 * 2\text{ml} = 210\text{ml} - 0,2\text{ml} = 210,2\text{ml}$$



$250\text{ml} - 1,5\text{div} = 250\text{ml} - 1,5 * 2\text{ml} =$
 $250\text{ml} - 3\text{ml} = 247\text{ml}$

$230\text{ml} - 0,4\text{div} = 230\text{ml} - 0,4 * 2\text{ml} =$
 $230\text{ml} - 0,8\text{ml} = 229,2\text{ml}$

$210\text{ml} + 0,1\text{div} = 210\text{ml} - 0,1 * 2\text{ml} =$
 $210\text{ml} - 0,2\text{ml} = 210,2\text{ml}$

$170\text{ml} + 3,4\text{div} = 170\text{ml} + 3,4 * 2\text{ml} =$
 $170\text{ml} + 6,8\text{ml} = 176,8\text{ml}$

oppure

$180\text{ml} - 1,6\text{div} = 180\text{ml} - 1,6 * 2\text{ml} =$
 $180\text{ml} - 3,2\text{ml} = 176,8\text{ml}$